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1. Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. NCCLS EP5-A Vol.19, No.2, February (1999).
2. Frolik CA. et al., Time-dependent changes in biochemical bone markers and serum cholesterol in ovariectomized rats: effects of raloxifene HCl, tamoxifen, estrogen, and alendronate. *Bone* (1996);18:621-7
3. Gallop PM. et al., Carboxylated calcium-binding proteins and vitamin K. *N Engl J Med*

(1980);302;1460-6.

4. Gaumet N. et al., Influence of ovariectomy and estradiol treatment on calcium homeostasis during aging in rats. *Arch Physiol Biochem* (1997);105:435-44.
5. Hauschka PV. et al., Direct identification of the calcium-binding amino acid, gamma-carboxy-glutamate, in mineralized tissue. *Proc Natl Acad Sci USA* (1975);72:3925-9.
6. Hauschka PV. et al., Osteocalcin and matrix Gla protein: vitamin K-dependent proteins in bone. *Physiol Rev* (1989);69:990-1047.
7. Kasugai Y. et al., Effects of tibolone (Org OD14) treatment for 3 months on ovariectomy-induced osteopenia in 8-month-old rats on a low-calcium diet: preventive testing for 3 months. *Bone* (1998);22:119-24. 8.
8. Lepola VT. et al., Bisphosphonates clodronate and etidronate in the prevention of ovariectomy-induced osteopenia in growing rats. *J Bone Miner Res* (1996);11:1508-17.
9. Poser JW. et al., Isolation and sequence of the vitamin K-dependent protein from human bone. Undercarboxylation of the first glutamic acid residue. *J Biol Chem* (1980);255:8685-91
10. Price PA. et al., Characterization of a gamma-carboxyglutamic acid-containing protein from bone. *Proc Natl Acad Sci USA* (1976);73:1447-51.
11. Rosenquist C. et al., Measurement of a more stable region of osteocalcin in serum by ELISA with two monoclonal antibodies. *Clin Chem* (1995);41:1439-45.
12. Seidlova-Wuttke D. et al., Pharmacology of Cimicifuga racemosa extract BNO 1055 in rats: bone, fat and uterus. *Maturitas* (2003);44Suppl.:S39-S50.

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生产者名称：英国 Immunodiagnostic Systems Limited  
生产者/生产场所地址：10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD, UK

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