# PRONTO® Usher / GSD™ kit

For t	he	dete	ectio	on of	fthe	follo	owir	١g	mut	atio	ns:

Usher Syndrome Type I..... R245X

Glycogen Storage Disease...... R83C

**REF** 9967

# Instructions for Use











#### INTENDED USE

The PRONTO<sup>®</sup> Usher / GSD kit is a Single Nucleotide Primer Extension ELISA Assay, intended for the qualitative *in vitro* detection of the following two mutations: R245X in the PCDH15 gene (Usher Syndrome Type 1) and R83C in the G6PC gene (Glycogen Storage Disease Type 1a), in amplified human DNA.

For in vitro diagnostic use.

### BACKGROUND

**Usher Syndrome** is an autosomal recessive disorder characterized by bilateral sensorineural deafness and progressive loss of vision due to retinitis pigmentosa. It is the most frequent cause of deafness and concurrent blindness with a prevalence of 1 in 16,000 to 1 in 50,000. There are three clinical subtypes of Usher syndrome, the most severe of which is Usher type 1 that involves deafness at birth, progressive blindness and balance problems. One of the Usher syndrome genes, PCDH15, is located on the long arm of chromosome 10. The R245X mutation accounts for a large proportion of cases of the type 1 Usher syndrome in Ashkenazi Jews, with a carrier rate as high as 1/40.

Glycogen storage disease type 1a (GSD1a, von Gierke disease) is an autosomal recessive disorder caused by a deficiency in the activity of glucose-6-phosphatase (G6Pase), a key enzyme in glucose homeostasis. Patients afflicted with GSD1a are unable to maintain glucose homeostasis and present with growth retardation, hypoglycemia, hepatomegaly, kidney enlargement, hyperlipidemia, hyperuricemia, and lactic acidemia. Long-term presentations include gout, hepatic adenomas with risk for malignancy, osteoporosis, platelet dysfunctionpulmonary hypertension, and renal failure.

GSD1a is caused by mutations in the G6-Pase gene, located on the long arm of chromosome 17. R83C is the most common mutation among Caucasians. It was found in 93% of Ashkenazi Jewish GSD1a patients, 71% of Mediterranean patients, 60% of Turkish patients, 48% of patients from South Europe and 28% of Hispanic patients. GSD1a has a disease rate of one in 100,000 births.

#### REFERENCES

- Lei et al, Genetic basis of glycogen storage disease type 1a: prevalent mutations at the glucose-6-phosphatase locus. Am. J. Hum. Genet. 57(4): 766-71 (1995).
- Parvari et al, Glycogen storage disease type 1a in Israel: biochemical, clinical, and mutational studies. Am. J. Med. Genet .72: 286-290 (1997).
- Ben-Yosef et al, A mutation of PCDH15 among Ashkenazi Jews with the type 1 Usher syndrome. N. Engl. J. Med. 48:1664-70 (2003).

## WARNINGS & PRECAUTIONS

- Reagents supplied in this kit may contain up to 0.1% sodium azide that is toxic if swallowed. Sodium azide has been reported to form explosive lead or copper azides in laboratory plumbing. To prevent the accumulation of these compounds, flush the sink and plumbing with large quantities of water.
- TMB Substrate solution is an irritant of the skin and mucous membranes. Avoid direct contact.
- The Stop Solution contains dilute sulfuric acid (1M), which is an irritant of the eyes and the skin. In case of contact with the eyes, immediately flush them with water. Do not add water to this product. In case of an accident or discomfort consult a physician (if possible, show the bottle label).
- In addition to reagents in this kit, the user may come in contact with other harmful chemicals that are not provided, such as ethidium bromide and EDTA. The appropriate manufacturers' Material Safety Data Sheets (MSDS) should be consulted prior to the use of these compounds.

### S ASSAY OVERVIEW

The PRONTO® procedure detects predefined polymorphisms in DNA sequences, using a single nucleotide primer extension ELISA assay.

- 1 TARGET DNA AMPLIFICATION: The DNA fragments that encompass the tested mutations are amplified. This amplified DNA is the substrate for the primer extension reaction.
- 2 POST-AMPLIFICATION TREATMENT: The amplified DNA is treated to inactivate free unincorporated nucleotides, so that they will not interfere with the primer extension reaction.
- **PRIMER EXTENSION REACTION:** A single-nucleotide primer extension reaction is carried out in a 96-well thermoplate. Each well contains a 5'-labeled primer that hybridizes to the tested DNA next to the suspected mutation site, and a single biotinylated nucleotide species, which complements the nucleotide base at the tested site. Each post-amplification treated sample is tested in two wells per mutation: the first well of each pair tests for the presence of the mutant allele (*mut*), while the second well tests for the

- presence of the normal allele (*wt*). The biotinylated nucleotide will be incorporated to the primer in the course of the reaction or not added, depending on the tested individual's genotype.
- 4 **DETECTION BY ELISA:** The detection of the biotinylated primers is carried out by an ELISA procedure: The biotin-labeled primers bind to a streptavidin-coated ELISA plate and are detected by a peroxidase-labeled antibody (HRP) directed to the 5' antigenic moiety of the primer. A peroxidase reaction then takes place in the presence of the TMB-Substrate.
- 5 INTERPRETATION OF THE RESULTS: Results are determined either visually (substrate remains clear or turns blue), or colorimetrically (substrate remains clear or turns yellow) following the addition of the stop solution.

#### O DISCLAIMER

- Results obtained using this kit should be confirmed by an alternative method.
- Confirmed results should be used and interpreted only in the context of the overall clinical picture. The manufacturer is not responsible for any clinical decisions that are taken.

The user of this kit should emphasize these points when reporting results to the diagnosing clinician or the genetic counselor.

MC9967 04.EN.02 Page 3 of 16 MC9967 04.EN.02 Page 4 of 16

### CONTENTS OF THE KIT

$\begin{array}{llllllllllllllllllllllllllllllllllll$
**
Solution D2 vials (red cap) (100 μL
ColoRed™ Oil1 dropper bottle (13 mL
Assay Solution
Wash Solution (conc. 20x)1 bottle (100 mL
HRP Conjugate1 vial (450 μL
TMB Substrate1 bottle (40 mL
Stop Solution (1M H <sub>2</sub> SO <sub>4</sub> )1 bottle (30 mL
PRONTO® Usher/GSD Plates*3 individually pouched plates
Detection Plates3 Streptavidin-coated ELISA plates

<sup>\*</sup> according to the customer's request

#### STORAGE AND STABILITY

- Store at 2-8°C. Do not freeze.
- Do not use the kit beyond its expiration date (marked on box label). Stability is maintained even when components are re-opened several times.
- Minimize the time reagents spend at room temperature.
- This kit has been calibrated and tested as a unit; do not mix reagents from kits with different lot numbers.

#### ADDITIONAL MATERIALS REQUIRED

- Tag DNA polymerase
- Deionized water (about two liters per kit)
- Thermowell plate or tubes (thin wall) for the post-amplification treatment
- Sterile pipette tips
- Troughs/reagent reservoirs for use with the detection reagents
- Thermocycler for a 96-well microplate
- Multichannel pipettes (5-50 µL and 50-200 µL)
- Positive displacement pipettes (1-5 μL, 5-50 μL, 50-200 μL & 200-1,000 μL)
- Filtered tips
- ELISA reader with 450 nm filter (optional 620 nm filter)
- Polaroid camera and color film to record results (optional)
- Automated microtiter plate washer, or a squirt bottle
- Vortex mixer
- Timer

# S ASSAY PROCEDURE

# 1 DNA AMPLIFICATION

- Dispense 2 μL template DNA (from an initial concentration of about 150 ng/μL) to a thermoplate well or tube.
- 2. **Prepare** a Master Mix in a sterile vial, according to the volumes indicated in the table below, plus one spare reaction volume. Add the Taq DNA polymerase to the amplification mix shortly before dispensing the mix. Mix gently by pipetting.

#### Master mix

Solution	Volume for one sample		
Amplification Mix Usher/GSD	13.0 µL		
Taq DNA polymerase (5 u/μL)	0.2 μL		

The following Taq DNA polymerases were validated for use with this procedure (lacking  $3' \rightarrow 5'$  exonuclease activity):

•	PHARMACIA	Cat. No. 27-0799
•	SIGMA	Cat. No. D-1806
•	ROCHE	Cat. No. 1-146-165
•	PROMEGA	Cat. No. M-1661
•	BIOLINE	Cat. No. M95801B
•	PERKIN ELMER	Cat. No. M801-0060
•	BIOLABS	Ca. No. M0267S

- 3. Dispense 13 µL Master Mix to each well or tube.
- Add one drop of ColoRed™ oil to each well. Do not touch the wells with the tip
  of the oil bottle. Even when using a thermocycler with a hot lid, it is essential to
  use oil.
- 5. **Place** the thermoplate well or tube in a thermocycler previously programmed with the following protocol:

### Cycling protocol

1.	94°C	5 min.	
2. 3. 4.	94°C 60°C 72°C	30 sec. 30 sec. 30 sec.	} 35 cycles
5.	72°C	5 min.	

MC9967 04.EN.02 Page 5 of 16 MC9967 04.EN.02 Page 6 of 16

 To verify amplification, subject 5 μL of the amplified product to electrophoresis in a 2% agarose gel.

#### Sizes of amplified fragments:

Gene	Mutation	Fragment Size
PCDH15	R245X	147 bp
G6PC	R83C	205 bp

#### Limitation of the test:

Different Taq DNA polymerases and thermocyclers may influence the amplification yield dramatically. Use a validated Taq DNA polymerase and a calibrated thermocycler.

# **2** POST-AMPLIFICATION TREATMENT

1

Only 10  $\mu\text{L}$  of each amplified DNA sample will be used to carry out this assay

1 **Prepare** a post-amplification treatment mix shortly before use. Combine in a single test tube the volumes appearing in the following table, multiplied by the number of tested samples, plus one spare volume.

#### **Post-Amplification Mix**

Solution	Volume for one sample
PRONTO® Buffer 2	30.0 μL
Solution C	1.0 µL
Solution D	0.75 μL

- **2 Mix** gently by pipetting this solution in and out five times. Do not vortex.
- 3 Add 30  $\mu L$  of the post-amplification mix into each well or tube containing 10  $\mu L$  of each amplified DNA sample.
  - Ensure that the solution you add becomes well mixed with the DNA sample by inserting the tip under the oil, down to the bottom of the tube and mixing the two solutions by pipetting.
- 4 Add one drop of ColoRed™ oil to each tube. Do not touch the tube with the tip of the oil bottle. Even when using a thermocycler with a hot lid, it is essential to use oil.
- **5 Incubate** for 30 minutes at 37°C, then for 10 minutes at 95°C in a thermocycler.

If not used immediately, the treated sample can be kept at 2-8°C for a maximum of four hours.

MC9967 04.EN.02 Page 7 of 16 MC9967 04.EN.02 Page 8 of 16

# 3 PRIMER EXTENSION REACTION

1 Program the thermocycler as follows:

Cycle		Temperature	Time	
Start:		94°C	15 sec.	
20 cycles:	ſ	94°C	30 sec.	
20 0yoloo.	ſ	60°C	10 sec.	
End:		18-25°C - Cool o	lown to roon	n temperature


**2 Take** a PRONTO<sup>®</sup> Plate out of its pouch. Notice the color at the bottom of the wells. For each mutation tested, use a pink well (*mut*) and a blue well (*wt*). Mark the plate with the ID numbers of your test.

If you intend to use less than a full plate, you can cut the plate and return the unused portion to the pouch. If you do this, seal the pouch immediately with its desiccant card inside.

3 **Dispense** 8 μL of post-amplification treated DNA into the first six wells in row A (see Fig. 1). Continue with the remaining samples. It is possible to transfer up to eight samples simultaneously using a multichannel pipette.

Ensure that the solution is at the bottom of each well by inspecting the plate from below. Be sure that the well does not contain air bubbles.

Figure 1: Scheme for dispensing Post-Amplification Treated DNA samples into the PRONTO<sup>®</sup> Usher/GSD™ Plate



#### Recommendation:

Use a new set of tips for each column. Alternatively use the same set of tips, but do not touch the bottom of the wells.

- **4** Tilt the plate and add one drop of ColoRed™ oil to each well. Do not touch the well with the tip of the oil bottle. Even when using a thermocycler with a hot lid, it is essential to use oil.
- **5** Turn on the thermocycler and start the cycling protocol.
- **6** When the thermal cycling is complete, you can proceed to the ELISA assay, or store the reaction products in the refrigerator and carry out the visualization steps within 24 hours.

# 4 ELISA ASSAY - COLOR DEVELOPMENT

The ELISA assay consists of the following steps:

- Binding the biotinylated primer to the streptavidin-coated plate.
- Washing away the unbound primer.
- Incubating with the HRP conjugate.
- Washing away the unbound conjugate.
- Incubating with the TMB substrate (color development).

The results of this assay can be determined in one of two ways:

**a** Visually: by monitoring the development of the blue color.

or

**Colorimetrically:** by adding Stop Solution and measuring the absorbance using an ELISA reader at a wavelength of 450 nm (yellow color).

Before proceeding with the ELISA assay, make your choice of visual or colorimetric determination of results.

# PREPARATION


- All components used in the detection step should reach room temperature before starting the assay.
- Dilute the 20x Wash Solution to 1x with deionized water.
   Dilute solution may be kept at 18-25° C for up to one month.
- **Peel** off the plastic cover of the Detection Plate. Mark the side of the plate with the kit name and test number.
- Place the PRONTO<sup>®</sup> Plate and the Detection Plate side by side, oriented in the same direction (see Fig. 2).

MC9967 04.EN.02 Page 9 of 16 MC9967 04.EN.02 Page 10 of 16

#### TRANSFER TO THE DETECTION PLATE

- 1 Fill a reagent reservoir / trough with the green colored Assay Solution. About 11 mL will be required for a 96-well plate.
- **2** Add 100  $\mu$ L of Assay Solution to the bottom of each well in column 1 of the PRONTO<sup>®</sup> Plate with a multichannel pipette. Mix the Assay Solution with the solution in the wells.
- 3 Without changing tips, transfer 100 μL from each well in this column to the first column in the Detection Plate (see Fig. 2).
  - Ensure that the solution at the bottom of all wells of the PRONTO<sup>®</sup> Plate has turned green by inspecting them from below.

**Figure 2:** Transferring the primer extension products from the PRONTO<sup>®</sup> Plate to Detection Plate.



- **4 Repeat** this procedure, using a new set of tips for each column. It is essential to maintain the order of the samples.
  - 10  $\mu$ L of oil carried over or 10  $\mu$ L of the sample left behind will not significantly affect the detection process.
- 5 Incubate for 10 minutes at room temperature (18-25°C).

#### DETECTION BY ELISA

		Visual	Colorimetric
Pro	ocedure	Detection	Detection
		(Blue color)	(Yellow color)
6	While the incubation with Assay Solution is		
	taking place, dilute the Conjugated HRP in	Dilution: 1:100	Dilution: 1:300
	Assay Solution .	(110 µL of	(37 µL of
	For every Detection Plate used (96 wells),	Conjugated HRP	Conjugated HRP
	about 11 mL of diluted conjugate is	into 11 mL	into 11 mL Assay
	required.	Assay Solution)	Solution)
	This solution should be freshly		
	prepared each time the test is run.		
7	Empty the plate and wash four times with		
	350 µL 1x Wash Solution. Ensure that the	$\sqrt{}$	$\sqrt{}$
	plate is relatively dry after the last wash	V	٧
	step.		
8	Add 100 µL of freshly diluted Conjugated		
	HRP to all the wells, with a multichannel	$\sqrt{}$	$\checkmark$
	pipette.		
9	Incubate at room temperature.	10 minutes	10 minutes
10	Wash the plate as in step 7.	V	√
11	Add 100 µL TMB-Substrate to each well		
	with a multichannel pipette and incubate at	15 minutes	15 minutes
	room temperature (18-25°C) until blue		
	color appears		
12	Add 100 µL of Stop Solution to each well		
	with a multichannel pipette. The solution	_	100 μL
	will turn yellow immediately.		
13	The results can be documented using a		
	Polaroid camera with color film (for	Agitate the plate	
	example - Fuji FP-100C), or by reading the	gently and read	_
	absorbance using an ELISA reader (signal	results at	
	wavelength setting).	O.D. 620 nm	
14	Within two hours read the absorbance		
	using an ELISA reader (single wavelength	_	O.D. 450 nm
1	setting)		

MC9967 04.EN.02 Page 11 of 16 MC9967 04.EN.02 Page 12 of 16

#### VALIDATION OF THE RESULTS


#### For Visual Detection:

For every mutation site tested, at least one of the wells should develop a deep blue color. Otherwise, the results are invalid for the relevant mutation (see Fig. 3).

#### For Colorimetric Detection:

For every mutation site tested, at least one of the two wells should yield an O.D. ≥0.50 reading.

Figure 3: Visual Interpretation of Genotypes




### INTERPRETATION OF RESULTS

**Important:** Heterozygous or homozygous mutant results should be confirmed by retesting. It is recommended to repeat the test with newly extracted DNA.

# **Criteria for Visual Interpretation**

A deep blue color indicates a positive signal, while negative signals appear as a clear or pale blue-colored well (see Fig. 4).

Figure 4: Genotype assignment according to visual inspection of test results.



# **Criteria for Colorimetric Interpretation**


The genotype of each sample is determined according to two criteria:

- 1. The O.D. values of the *mut* and *wt* wells.
- 2. The ratio of mut / wt O.D. values.

Calculate the *mut / wt* ratios by dividing the signal of the *mut* well by the signal of the *wt* well.

Identify the correct genotype using the table below:

Genotype	mut well	wt well	mut/wt ratio
	(O.D. 450)	(O.D. 450)	
Normal	O.D. <u>&lt;</u> 0.35	O.D. <u>&gt;</u> 0.5	ratio < 0.5
Heterozygote	O.D. <u>≥</u> 0.5	O.D. <u>&gt;</u> 0.5	0.5 < ratio < 2.0
Homozygote	O.D. <u>&gt;</u> 0.5	O.D. <u>&lt;</u> 0.35	ratio > 2.0



MC9967 04.EN.02 Page 13 of 16 MC9967 04.EN.02 Page 14 of 16

# Usher/GSD® - PROCEDURE SUMMARY

DNA EXTRACTION: from human whole blood, using a validated method.

#### DNA AMPLIFICATION:

Volumes per reaction: 2 µL Template DNA + 13 µL Amplification Mix + 0.2 µL Taq Polymerase.

Cycling protocol:

94°C 5 min→35 cycles of {94°C 30 sec. / 60°C 30 sec. / 72°C 30 sec.} →72°C 5 min.

#### POST-AMPLIFICATION PROCEDURE:

PRONTO<sup>®</sup> Buffer 2 Volumes per reaction: 30.0 µL

Solution C 1.0 µL Solution D  $0.75 \, \mu L$ 

■ Pipette in and out to mix.

■ Add 30 µL into each well containing 10 µL amplified product, mix well.

■ Add one drop of ColoRed<sup>TM</sup> oil.

■ Incubate 30 minutes at 37°C, then 10 minutes at 95°C.

#### PRIMER EXTENSION REACTION:

Dispense 8 µL of each post-amplification treated DNA into six wells of the PRONTO<sup>®</sup> Plate.

■ Add one drop of ColoRed<sup>TM</sup> oil.

Start the cycling protocol:

94°C 15 sec→20 cycles of {94°C 30 sec. / 60°C 10 sec.} →Cool.

#### **DETECTION:**

 Add 100 µL Assay Solution to each well in the PRONTO<sup>®</sup> Plate and mix.
 Transfer 100 µL from each well of the PRONTO<sup>®</sup> Plate to the identical position in the Detection Plate. Incubate 10 minutes at RT.

**Empty** the wells and wash four times with 350 uL of 1x Wash Solution.

For either visual or colorimetric detections, continue as follows:

	Visual Detection	Colorimetric Detection
Add 100 µL of Conjugated HRP to every well and incubate for 10 minutes at RT.	Dilution 1:100	Dilution 1:300
Empty the wells and wash four times with 350 µL of 1x Wash Solution.	√	<b>V</b>
Add 100 µL of TMB Substrate to each well and incubate at RT for:	15 minutes	15 minutes
Add Stop Solution	_	100 μL per well
Read O.D. at:	620 nm	450 nm

For troubleshooting guide, please refer to our website: www.prontodiagnostics.com/ts

PRONTO is a Registered Trademark COLORED is a trademark of Pronto Diagnostics Ltd. Copyright © Pronto Diagnostics Ltd.

The Pronto<sup>®</sup> Technology is covered by US patent 5.710.028, by European patent 0648222 and by corresponding national patents.

### Produced by Pronto Diagnostics Ltd.

Kirvat Atidim Building 3, 3<sup>rd</sup> floor POB 58129 Tel Aviv 6158002 Israel

Tel: +972.73.2126155 Fax: +972.73.2126144

**Customer Service:** E-Mail: info@prontodiagnostics.com

#### **Authorized EU Representative:**

MedNet GmbH. Borkstrasse 10, 48163 Münster, Germany

Tel: +49.251.32266-0 Fax: +49.251.32266-22

MC9967 04.EN.02

MC9967 04.EN.02 MC9967 04.EN.02 Page 15 of 16 Page 16 of 16