



TECHNICAL SHEET

PRODUCT NAME: Anti-S100A12 pAb (protein A purified)

Code:	MT-20HSM
Clone:	NA
Specie:	Rabbit pAb
Physical State:	Liquid
Host Animal (source):	Rabbit pAb
Specificity:	20HSM is specific for S100A12 detection.
Cross Reactivity:	NA
Purification:	Anti-S100A12 pAb has been purified in a multi-step procedure including a protein A chromatography $\geq 95\%$ pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is lower than 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 14/01/2019