



TECHNICAL SHEET

PRODUCT NAME: Anti-S100b pAb (Affinity purified)

Code: MT-20HSBA

Clone: NA

Physical State: Liquid

Host Animal (source): Rabbit pAb

Specificity: 20HSBA is specific for S100b detection.

Cross Reactivity: NA

Purification: Anti-S100b pAb has been purified in a multi-step procedure. $\geq 95\%$ pure by SDS-PAGE.

Storage Conditions: Storage at 2-8°C.

Health & Safety Information: This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

Presentation: PBS pH: 7.4, 0.09% Sodium Azide (NaN_3).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 19/07/2019