

TECHNICAL SHEET

PRODUCT NAME: SC2 recombinant Spike Glycoprotein (S) Brazilian variant

Code: MT-25C19S.p.1, MT-25C19S.p.1s1

Physical State: Liquid

Source: Hek293 cells

Description: SC2 recombinant Spike Glycoprotein (S) Brazilian variant

expressed and purified from Hek293 cells.

Appearance: Clear

Purity: >90%

Storage Conditions: Short term storage at 2-8°C. Long term storage at -20°C. Avoid

multiple/freeze thaw cycles by storing multiple aliquots at

-20ºC.

Shipping: With cool packs.

Health & Safety Information: Good Laboratory Practices should be followed when handling

this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in

accordance with the respective national regulations.

CerTest

BIOTEC

Presentation: Tris buffer with additives and salts, pH: 8.00.

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 28/04/2021

